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Intensiv HygienicTray

Product description

- Easily dismantled in 2 parts
- Suitable for the standard sterilisation tray
- Ideal arrangement of the 13 FG and 5 RA perforations for easy access to the instruments
- Stainless steel
- Dimensions
H x L x W: 30 x 49 x 29mm

Indications

- HygienicTray may be used as an organizer for a specific indication treatment, i.e. prep bur block. In addition, it enables an easier handling of the instruments during the treatment and provides a support during the sterilisation process.

Practical recommendations

- To avoid any cross contamination, use disinfected tweezers to take the instruments.

Maintenance and sterilisation

- For easy cleaning, dismantle the Intensiv HygienicTray completely (cover, instrument tray).
- The HygienicTray is supplied non-sterile. Prior to first use and immediately after each use, the HygienicTray must be disinfected, cleaned, and sterilized.
- During the sterilisation keep the HygienicTray separate from non-stainless steel instruments, such as rubber polishers and abrasive stones.
- If disinfection by immersion is being used, it is important to use only cleaning/disinfection solutions that provide corrosion protection and strictly observe the concentrations recommended by the manufacturer. Sterilization must be performed with validated procedures. If possible, use a single-pulsed or fractionated vacuum autoclave and subvacuum drying. Chemiclave sterilizers may also be used.

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